
Pedal Exerciser (model 740-742)

Contents:

- 2 Pedal mounts
- 1 Pedal assembly
- 4 Allen head bolts
- 4 7/16" self locking nuts
- 1 Rubber gasket
- 1 Tension knob
- 1 Allen key

Assembly Instructions:

1. Line up the pedal mounts so they are facing each other.
2. Place the rubber gasket on top of the pedal mounts, lining up the holes for the bolts.
3. Place the pedal assembly on top of the rubber gasket, again lining up the holes for the bolts.
4. Insert the Allen head bolts through the pedal assembly, the rubber gasket and the pedal mounts and attach the nuts (nylon end facing down) from underneath.
5. Check to see if the pedal mounts are even. Then, using a 7/16" wrench or pliers to hold the nuts in place, tighten the hex bolts securely with the Allen key.
6. Install the tension knob as illustrated (it controls the tension of the Pedal exerciser).

Usage Instructions:

Enjoy a bicycle-like exercise by using your ProActive™ Pedal Exerciser while sitting in your favorite chair, on the side of a bed or in a wheelchair. It can also be placed on a table in front of you to exercise your arms and upper body.

You can easily increase the level of resistance (and of desired effort) by turning the tension knob clockwise. Turn it counter-clockwise to ease the tension and make pedaling easier.

Start slowly and at a comfortable tension level, gradually building up to a goal of 20 min 3 times per week. As you progress, increase the tension to continue getting maximum benefits from your workout. Do not attempt to do too much right away as it may result in injury and/or soreness.

IMPORTANT

Check the tightness of the bolts periodically. The bolts may require tightening from time to time. Do not stand up on the pedals. Use only while seated. Always consult a physician before starting a new exercise program.

Questions ? Comments ?
1-800-363-2381

AMG Medical Inc.
Montreal, QC H4T 1V5
West Chazy, NY 12992
www.amgmedical.com

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