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Montréal, QC H4T 1V5 1-800-363-2381 West Chazy, NY 12992 1-888-412-4992

[www.proactive-painmanagement.com](http://www.proactive-painmanagement.com)

## INSTRUCTIONS: SHOULDER AND NECK COMPRESS

Compress should be applied for 20 minutes, or as recommended by a doctor. Remove and wait at least 20 minutes. Repeat as needed. To further reduce hot or cold level, place compress in soft cloth covering before applying.

### For Heat Therapy:

Microwave compress per the below suggested heating times and test compress temperature. Heat for additional 5 seconds as needed until desired temperature is reached.

#### Suggested heating time for one (1) compress:

700 Watts 60 seconds	1000 Watts 50 seconds	1250 Watts 40 seconds
-------------------------	--------------------------	--------------------------

#### Do not exceed heating time for one (1) compress:

700 Watts 80 seconds max.	1000 watts 70 seconds max.	1250 Watts 60 seconds max.
------------------------------	-------------------------------	-------------------------------

### For Cold Therapy:

Freeze for at least 2 hours in a freezer storage bag.

### Caution:

Always test compress temperature prior to use. Overheating can cause damage or destroy it. Prolonged use may cause injury to user including burns. For external use only. Do not use on infants. Never sit, lie, or lean against the compress. Do not puncture compress or ingest contents. Do not apply to sensitive/insensitive skin or open wounds. Use only for intended purposes. Do not fall asleep using compress. Use only under adult supervision. For chronic, extreme or prolonged pain, please contact your physician.

### Storage:

Always keep compress in a freezer storage bag when kept in the freezer.

### Cleaning:

To clean compress, use mild soap and water and air dry before reusing.

Manufactured and printed in China for:

A.M.G. Medical Inc.

Montréal, QC H4T 1V5 1-800-363-2381 West Chazy, NY 12992 1-888-412-4992

[www.proactive-painmanagement.com](http://www.proactive-painmanagement.com)

Rev 1 0417 P/N 911-740808