



**NEW!**

# Elastitube™

Now, the brand professionals trust offers you a choice line of tubular stockinette bandages.



### Skin Friendly

- Latex-free to help prevent skin irritation.
- Remains secure on wounds and sensitive areas to help protect skin.

### Autoclavable

- May be processed to provide sterile dressings.

### Water-Repellent

- Woven, continuous tubular fabric made of water-repellent, top quality cotton yarn.

### Practical

- Preferred wound care alternative to using flat gauze bandages, since no winding is required to be secured in place and to permit freedom of movement.



### Featuring user-friendly packaging

- Size indicated on multiple surfaces for easy identification.
- Push tab is easy to perforate with fingers to create dispensing slot.
- Dispensing side features a body silhouette with colored areas that represent where Elastitube size can be applied.
- A summary table of all sizes, along with their corresponding body part coverage, appears on the back panel.



**NEW!**

# Elastitube™

**Ideal for a wide range of applications:**

- Acts as skin shield by protecting skin beneath plaster casts or plaster of Paris bandages.
- Used in general and plastic surgery as well as dermatology.
- Suitable as a traction bandage to cover splints and cotton wool padding.

Size	Item No.	Body Parts	Approx. Length	Approx. Width*
1	118-571	Fingers and toes	21.9 yd / 20 m	0.59" / 1.5 cm
2	118-572	Large fingers and toes	21.9 yd / 20 m	0.98" / 2.5 cm
3/4	118-573	Child's arm, finger <b>with splints (adult)</b>	21.9 yd / 20 m	1.38" / 3.5 cm
5/6	118-575	Child's leg, <b>adult</b> arms and thin legs	21.9 yd / 20 m	2.36" / 6 cm
7/8	118-576	Child's head, <b>medium adult</b> legs	21.9 yd / 20 m	3.15" / 8 cm
9	118-577	<b>Adult</b> head and <b>large</b> legs	21.9 yd / 20 m	3.94" / 10 cm

\* Unstretched. Actual width may vary slightly due to the effect of environmental factors on the product.

**Elastitube is made of 100% pure cotton and is latex-free.**



To learn more about MedPro® Elastitube™ or any other AMG product, contact your AMG Representative, call toll-free at 1-800-363-2381 or visit [www.amgmedical.com](http://www.amgmedical.com).

**AMG** AMG Medical Inc.  
Montreal, QC H4T 1V5  
[www.amgmedical.com](http://www.amgmedical.com)



REV. 1.0 1104