

Proactive - Product Comparison Sheet



TENS	Thera3 715-430	TheraHeat 715-440	Alevia 715-425	Pulse 715-420	Libera 715-435
TENS	Yes	Yes	Yes	Yes	Yes
EMS	Yes	No	Yes	No	No
Massage	Yes	No	No	No	No
Heat	No	Yes	No	No	No
Wireless	No	Yes	No	No	Yes
Presets	78	8 pulse/2 heat	18	3	5
Channels	Dual	Single	Dual	Dual	Single
Customized Programming	Yes	No	Yes	Analog Manual	No
Rechargeable	No	Yes USB Cable Included	No	No	Yes USB Cable Included
Power Supply	4 X AA Batteries	3.7V Battery	9V Battery Included	9V Battery Included	3.7V Battery
Warranty	1 Year Limited	1 Year Limited	1 Year Limited	1 Year Limited	1 Year Limited